

Audéo ZIP

Product information



Audéo ZIP features state-of-the-art CORE performance features in a tiny, InstantFit CIC housing. Fitted with a wide range of domes, Audéo ZIP can be configured as a discreet fully automatic system or, by teaming it with a KeyPilot, provide access to additional manual programs to suit a busy

lifestyle. The ultimate in cosmetics, Audéo ZIP is specifically for first time users who are looking for an in-the-ear solution which they can try on the spot.

Audéo ZIP is available in three performance/price levels.

Features/Functions

	Audéo ZIP IX	Audéo ZIP V	Audéo ZIP III
SoundFlow	Premium	Advanced	Standard
Manual programs	3	3	2
WhistleBlock Technology	■	■	■
SoundRecover	■	■	■
Directionality	omnidirectional	omnidirectional	omnidirectional
SoundRelax	■		
EchoBlock System	■		
WindBlock Management	■	■	
NoiseBlock Processing	20 channels	16 channels	6 channels
DataLogging	■	■	■

Fitting software


- iPGF version 2.6 or higher, NOAH compatible

Hardware


- PC (IBM compatible) with NOAHlink and/or HI-PRO interface

Programming

Connecting the programming cable to an Audéo ZIP



1. Open the battery compartment and lift the cover out of the hinge, using the battery compartment tool (twisting motion).






2. Connect the cable. Attach the cables with the blue and red markers to the left and the right hearing instruments, respectively.




PHONAK life is on

Key features





The complete range of data sheets are available in www.phonakpro.com

		Audéo ZIP IX	Audéo ZIP V	Audéo ZIP III
				
Performance profile				
Maximum Power Output (dB SPL)	2 cc coupler	112	112	112
	Ear Simulator	119	119	119
Maximum gain (dB)	2 cc coupler	40	40	40
	Ear Simulator	50	50	50
Frequency range - Ear Simulator (Hz)		150 Hz - 7900 Hz	150 Hz - 7900 Hz	150 Hz - 7900 Hz
Working current (mA)		1.0 mA	1.0 mA	1.0 mA
Battery size		10	10	10

Fitting range



Options

	Audéo ZIP IX	Audéo ZIP V	Audéo ZIP III
 Open dome (SDS 2.0; 3 sizes: small, medium, large)	S	S	S
 Closed dome (SDS 2.0; 3 sizes: small, medium, large)	0	0	0
 Power dome (SDS 2.0; 3 sizes: small, medium, large)	0	0	0
 KeyPilot	0	0	0
Choice of 3 colors: Cocoa/Black transparent (H1), Pink/Black transparent (L5), Brown/Black transparent (L6)	0	0	0

Note:

Some options are country specific.
Please refer to your local order form.

S = Standard

0 = Optional

SDS = Sound Delivery System

Technical description

- 1 Dome
- 2 Receiver housing with left/right indicator
- 3 MultiFlex Canal
- 4 SecureFit Tab
- 5 Microphone with protector
- 6 Removal pin
- 7 Battery compartment with on/off switch
- 8 Serial number (inside battery compartment)

